

# MICROSUTURE

MICROSUTURE is a nylon monofilament nonabsorbable, sterile surgical suture that is indicated for use in general soft tissue approximation and/or ligation, including use in microsurgical, cardiovascular, ophthalmic, and neurological procedures.

## Features & Benefits

### Characteristics

- Uniform diameter of the suture
- Maximum flexibility of the strand
- Consistent knot security
- Minimal memory



## Ordering Information

Code	Size	Needle Type	Chord	Diameter	Curvature	Needle Length (arch)	Suture Length	Picture
GEM192BK	8-0 Nylon	taper	5 mm cord	140 dia	3/8 curve	6.4 mm	15 cm/6" suture lgth	
GEM395BK	8-0 Nylon	tapersharp	5 mm cord	140 dia	3/8 curve	6.4 mm	15 cm/6" suture lgth	
GEM177BK	9-0 Nylon	taper	3 mm cord	100 dia	3/8 curve	3.8 mm	15 cm/6" suture lgth	
GEM180BK	9-0 Nylon	taper	4 mm cord	100 dia	3/8 curve	5.1 mm	15 cm/6" suture lgth	
GEM386BK	9-0 Nylon	tapersharp	3 mm cord	100 dia	3/8 curve	3.6 mm	15 cm/6" suture lgth	
GEM389BK	9-0 Nylon	tapersharp	4 mm cord	100 dia	3/8 curve	5.1 mm	15 cm/6" suture lgth	
GEM171BK	10-0 Nylon	taper	3 mm cord	70 dia	3/8 curve	3.8 mm	10 cm/4" suture lgth	
GEM385BK	10-0 Nylon	tapersharp	3 mm cord	100 dia	3/8 curve	3.8 mm	10 cm/4" suture lgth	
GEM170BK	11-0 Nylon	taper	3 mm cord	70 dia	3/8 curve	3.8 mm	10 cm/4" suture lgth	
GEM195BK	11-0 Nylon	taper	2.5 mm cord	50 dia	3/8 curve	3.6 mm	10 cm/4" suture lgth	
GEM382BK	11-0 Nylon	tapersharp	3 mm cord	70 dia	3/8 curve	3.8 mm	10 cm/4" suture lgth	

### Black Needle

GEM192BKC	8-0 Nylon	taper	6 mm cord	140 dia	3/8 curve	6.4 mm	15 cm/6" suture lgth	
GEM395BKC	8-0 Nylon	tapersharp	6 mm cord	140 dia	3/8 curve	6.4 mm	15 cm/6" suture lgth	
GEM180BKC	9-0 Nylon	taper	5 mm cord	100 dia	3/8 curve	5.1 mm	15 cm/6" suture lgth	
GEM389BKC	9-0 Nylon	tapersharp	5 mm cord	100 dia	3/8 curve	5.1 mm	15 cm/6" suture lgth	

### INDICATIONS FOR USE:

MICROSUTURE is indicated for use in general soft tissue approximation and/or ligation, including use in microsurgical, cardiovascular, ophthalmic, and neurological procedures.

**Rx Only.** For safe and proper use of this device, refer to the Instructions for Use.

Manufactured by Riverpoint Medical and Distributed by Synovis Micro Companies Alliance, Inc.